**Instructions for use**

**Ambu<sup>®</sup> PEEP 10 valve and  
Ambu<sup>®</sup> PEEP 20 valve**

<b>English</b>	
Instructions for use .....	1
<b>България</b>	
Инструкции за употреба .....	9
<b>Česky</b>	
Návod k použití .....	17
<b>Dansk</b>	
Brugervejledning .....	25
<b>Deutsch</b>	
Bedienungsanleitung .....	33
<b>Ελληνικά</b>	
Οδηγίες χρήσης .....	41
<b>Español</b>	
Instrucciones de Uso .....	49
<b>Eesti</b>	
Kasutusjuhend .....	57
<b>Suomi</b>	
Käyttöohjeet .....	65
<b>Français</b>	
Mode d'emploi .....	73
<b>Hrvatski</b>	
Upute za upotrebu .....	81
<b>Magyar</b>	
Használati útmutató .....	89
<b>Italiano</b>	
Istruzioni per l'uso .....	97
<b>日本語</b>	
使用方法 .....	105
<b>Lietuviškai</b>	
Naudojimo instrukcija .....	113
<b>Latviski</b>	
Lietošanas pamācība .....	121
<b>Nederlands</b>	
Gebruiksaanwijzing .....	129
<b>Norsk</b>	
Bruksanvisning .....	137
<b>Polski</b>	
Instrukcja użycia .....	145
<b>Português</b>	
Instruções de utilização .....	153
<b>Română</b>	
Instrucțiuni de utilizare .....	161
<b>Русский</b>	
Инструкция по применению .....	169
<b>Slovensky</b>	
Návod na použitie .....	177
<b>Slovenščina</b>	
Navodila za uporabo .....	185
<b>Svenska</b>	
Bruksanvisning .....	193
<b>Türkçe</b>	
Kullanım talimatları .....	201
<b>中文</b>	
使用说明 .....	209

<b>Contents</b>	<b>Page</b>
1 Warnings/Precautions	2
2 Introduction	3
3 Specifications	3
4 Description	4
5 Operating instruction	5
6 Cleaning – disinfecting – sterilisation	5
7 Trouble–shooting	6
8 Spare parts	7

**Note:** These directions for use may be updated without further notice.  
Copies of the current version are available from the manufacturer.

Ambu® are registered trademarks of Ambu A/S, Denmark.  
Ambu A/S is certified according to ISO 9001 and ISO 13485.

## Directions for use

# Ambu<sup>®</sup> PEEP 10 valve and Ambu<sup>®</sup> PEEP 20 valve

## 1. Warnings/Precautions

### WARNING

- The PEEP valve should only be used by trained persons who are aware that PEEP may have an adverse effect on the patient's circulation.
- The PEEP 20 valve should only be used in situations where the patient's condition can be checked carefully by means of bloodgas and haemodynamic tests.
- Do not use the PEEP valve on ventilators or breathing systems if the attachment of such valves is not specified by the manufacturer.

### CAUTION

- U.S. federal law restricts this device to sale by or on the order of a physician (USA and Canada only).
- Following disassembly and cleaning of the PEEP valve it is important to check that the parts are in perfect condition, e.g. that there are no cracks, and that the internal valve assembly does not leak silicone oil, and to assemble them correctly.
- For the adjustment of the PEEP 20 valve, a manometer should be connected to the patient breathing for monitoring of the adjustment. For the PEEP 10 valve, the same applies when using the PEEP 10 valve on infants (below 10 kg body weight).

Increase PEEP during ventilation  
to min. required level

## 2. Introduction

The Ambu® **PEEP 10 valve** and Ambu® **PEEP 20 valve** are suitable for use with:

- Resuscitators
- Ventilators, if specified by the manufacturer.
- CPAP systems, if specified by the manufacturer.
- When used on resuscitators the PEEP valve influences neither the inspiratory O<sub>2</sub> concentration nor the inspiratory resistance of the system. The PEEP valve can therefore remain at the patient valve during both spontaneous breathing and resuscitation.


## 3. Specifications

The Ambu® **PEEP 10** and **20** valves are in conformity with Council Directive 93/42/EEC concerning Medical Devices.

### Adjustment range:

PEEP 10: 0.15-1.0 kPa (1.5-10 cm H<sub>2</sub>O)

PEEP 20: 0.15-2.0 kPa (1.5-20 cm H<sub>2</sub>O)

*Characteristic curves of the pressure drop through the valves.*

Dimensions: Length 65 mm, diameter 45 mm

Weight: 40 g

Recommended temperatures:


Operation: -20 °C to +50 °C (-4 °F to +122 °F)

Storage: -40 °C to +70 °C (-40 °F to +158 °F)

<b>Parts</b>	<b>Materials</b>
Transparent plastic parts	Polysulphone
Plastic parts of valve body	Polysulphone with delrin disc and silicone oil
Adjustment cap	Polyamide (nylon)
Sealing in adjustment cap	EPDM rubber
Springs	Stainless steel

## 4. Description

The PEEP valve consists of three main parts:

## 5. Operating instruction

After connecting the PEEP valve to the equipment adjust the valve to the required PEEP level by turning the cap.

**Note:** The cap is not designed to be unscrewed from the housing. When the cap is turned counterclockwise towards the zero setting resistance is felt to prevent the cap from being unscrewed by accident. It does, however, not prevent unscrewing of the cap and once this is attempted, the locking ability is weakened.

For the adjustment of the PEEP 20 valve, a manometer should be connected to the patient breathing for monitoring of the adjustment. For the PEEP 10 valve, the same applies when using the PEEP 10 valve on infants (below 10 kg body weight).

### 5.1. Disposal of product

Used products must be disposed of according to local procedures.

## 6. Cleaning – disinfecting – sterilisation

The PEEP valve is exposed to expiratory gases. It should, according to local disinfection instructions, be cleaned and disinfected after each patient. Before cleaning, disassemble the three parts shown in section 4. No further disassembly is necessary.

### 6.1. Cleaning of parts

Wash the parts in warm, soapy water, rinse thoroughly in tap water and let the parts dry before reassembly.

An automatic washing machine for anaesthesia equipment may be used.

## 6.2. Disinfecting

Chemical disinfectant can be used. The supplier's directions should be followed carefully and the disinfectant must be suitable for the materials (paragraph 3).

## 6.3. Sterilisation

The parts can be autoclaved at 134 °C (273 °F) after cleaning.

## Reassembly

Ensure that the ends of the central spindle of the internal valve assembly are mounted correctly both in the hole of the inlet connector and in the housing with the adjustment cap.

# 7. Trouble-shooting

**Problem:** It is not possible to obtain correct adjustment of the required PEEP level.

**Possible cause:** Leak in the system.

**Correction:** Check all the connections systematically.

**Problem:** Valve noise.

**Possible cause:** Insufficient damping.

**Correction:** Disassemble the valve, move the valve disc back and forth in relation to the valve spindle. The valve should be serviced if this does not help.

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